

# 4PT™

4 POLE TECHNOLOGY



## Interceptor S2

  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

# Interceptor S2



## Features

- Brand new sporty streamline design
- Fully adjustable suspension system that delivers a smooth and comfortable ride
- Modern LED lighting system delivers a brilliant light with less power consumption
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Built-in half speed switch which operates automatically when turning to provide extra safety
- With Merits 4PT (4 pole technology) and large controller delivering much more power for climbing steep driveways and hills
- With new easy to access 12 volt socket



Modern LED lighting system delivers brilliant light with less power consumption




More space for storage

Digital LCD dashboard display

Handy 12 volt socket

## Dimensions cm/in



Overall Dimension : A x B x C 154 x 72 x 145.5cm/61 x 29 x 57.3"

Seat Dimension : D x E x F 51 x 46 x 63cm/20 x 18 x 25"

Seat Height from Ground : G 67.1cm/26.4"

Seat Height from Deck : H 48.6cm/19"

Wheel Base : I 113cm/45"

## Weight kg/lb

Base	Battery/ea	Seat	Total with Battery	Capacity
84/185	25/55	18/40	165/364	204/450

## Technical Specs

Max Speed	14/8.9mph
Range	40-60km/26-38mi
Turn Radius	180cm/70.9"
Ground clearance	13cm/5.1"
Motor	DC 24V, 650W (Rating)/2200W (Max)
Controller	PG S-Drive 200A
Battery	12V/60-100A/H x 2pcs
Charger	8A off-board
Gradient	10°
Front Wheel	14" (3.50-8) pneumatic
Rear Wheel	16" (4.00-10) pneumatic
Brake	Intelligent, regenerative, electromagnetic brakes

Australian Standards Compliant AS/NZS 3695.2:2013

  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.