



REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION
1	64210011	CONTROLLER ONLY, SHARK-50A, DK-PMB50	14	32600498	REAR WHEEL ASS'Y., 3.00-8 (4PR/C248 GREY)	27	32222085	CAPTAIN SEAT ONLY, WITH/MC7/15/22/25 COLOR 2/TONE 20 SW*18D*20H W/BE
2	43700209	VELCRO TAPE (M), JET 1	15	42500110	ALUM 3.00-8 RIM	28	32100180	SEAT MOUNTING ASS'Y
3	43700210	VELCRO TAPE, FEMALE 50.8W*300L	16	42500111	ALUM 3.00-8 IMMER RIM	29	41300416	EXTRUSION, ALUMNUM, HOOK FOR SEAT
4	43600540	CLIP, SNAP-IN POM φ6.3*25.4L	17	42500112	HUB φ17*36.5L	30	31000078	ARMREST ASS'Y., LH
5	54200011	HITCH-PIN 8MM OD, 40MML	18	43600011	3.00-8 (4PR/C248 GREY) FOAM FILLED TIRE	31	31000088	ARMREST ASS'Y., RH
6	44600508	SHROUD ONLY, BLACK, LH., P327	19	30510899	FENDER, CANDY APPLE RED, LH	32	43700319	SEAT BELT 1725MM LOOP TYPE
7	44600509	SHROUD ONLY, BLACK, RH., P327	20	30510900	FENDER, CANDY APPLE RED, RH	33	33503064	JOYSTICK BRACKET ASS'Y., SHARK, ED
8	42500031	ALUM. FOOTPLATE, BLACK	21	30500173	SHROUD ASS'Y.	34	64810013	JOYSTICK ONLY, DY SHARK, DK-REMD01
9	33100728	CASTOR FORK ASSY, HD(50mm),8"*2" PU WHEEL, GRAY	22	43600568	PLASTIC FASTEN BOLT, BLACK,1/4",DIA 6.3*21.6L	35	63300042	ELECTRONIC, HARNESS, DY FLIGHT BUS CABLE, 1.5M,
10	33000130	WHEEL ASS'Y., 200*50 PU, GRAY.	23	43600608	REAR SHROUD ONLY, BLACK	36	65000084	CHARGER
11	41700223	ABSORBER 165L 300LBS/IN	24	32102001	TRAPEZE BAR UPRIGHT WELDMNT, 2 HOLE, FRONT	37	63200088	BATTERY CABLE
12	62012021	MOTOR & GEARBOX ASS'Y.,LH., M9L01EA-5	25	32103073	SEAT TOWER, HORIZONTAL, FRONT TRAPEZE BAR, ED.	38	63334601	WIRE, CONTROLLER, 440L, P326
13	62022021	MOTOR & GEARBOX ASS'Y.,RH., M9L01EA-5	26	32100076	SEAT TOWER ASS'Y, HORIZONTAL, REAR TRAPEZE BAR			