



MODEL#S940A

Rev: 0



REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION
1	30300338	FRONT SUSPENSION ASSY(L)	23	30541112	FRONT SHROUD ASSY, BRIGHT RED	35	30561038	FRONT WINDSHIELD FOLLOWING COVER ONLY, ORANGE
2	30300339	FRONT SUSPENSION ASSY(R)	24	44600591	RUBBER MAT, BLACK	35	30550051	FRONT WINDSHIELD FOLLOWING COVER ONLY, SILVER DEEP GREEN
3	41700290	SUSPENSION(HYDRAULIC),165L,150LBS/IN, YELLOW	25	30500159	STORAGE COVER	35	30541080	FRONT WINDSHIELD FOLLOWING COVER, BRIGHT RED
4	30303101	FRONT WHEEL AXLE ONLY, LH	26	30519489	SHROUD ASSY, LH, SILVER WHITE	36	43690374	FRONT WINDSHIELD REAR COVER
5	30303102	FRONT WHEEL AXLE ONLY, RH	26	30561076	SHROUD ASSY, LH, ORANGE	37	43601374	FRONT WINDSHIELD COVER
6	30302009	ADJUST SHAFT, LONG, ASSY	26	30550111	SHROUD ASSY, LH, SILVER DEEP GREEN	38	43690379	RUBBER MAT, BLACK
7	30302017	ADJUST SHAFT, SHORT φ10*175L(ZP)	26	30541084	SHROUD ASSY, LH, BRIGHT RED	39	61200174	HEAD LIGHT ASSY(L)
8	43690177	FRONT BUMPER, S940	27	30519490	SHROUD ASSY, RH, SILVER	40	61200175	HEAD LIGHT ASSY(R)
9	32100025	UPPER SEAT TUBE ASSY.	27	30561075	SHROUD ASSY, RH, ORANGE	41	30400391	REAR FRAME ASSY,
10	61120085	CIRCUIT BREAKER, 100A	27	30550112	SHROUD ASSY, RH, SILVER DEEP GREEN	42	41700291	SUSPENSION(HYDRAULIC),210L,150LBS/IN, RED
11	63745602	BREAKER WIRE	27	30541085	SHROUD ASSY, RH, BRIGHT RED	43	61090063	SPRING SWITCH, PS8-23
12	63748101	HARNESS, HALF SPEED SWITCH	28	30519558	REAR LAMPS OF TOP COVER ONLY, SILVER WHITE	44	62009008	TRANSAXLE & MOTOR ASSY., M817AG-21
13	30303336	TURN TUBE COMPLETE (ED)	28	30561409	REAR LAMPS OF TOP COVER ONLY, ORANGE	45	62203026	TRANSAXLE (M817)
14	33100724	FRONT WHEEL ASSY, 14*3.50-8, BLACK	28	30550135	REAR LAMPS OF TOP COVER ONLY, SILVER DEEP GREEN	46	62100160	MOTOR ASSY
15	42500099	INNER RIM ONLY, FOR 14*3.00-8	28	30541113	REAR LAMPS OF TOP COVER ONLY, BRIGHT RED	47	33200153	REAR DISK BRAKE ASSY
16	42500257	OUTTER RIM ONLY, FOR 14*3.00-8 FRONT WHEEL	29	30500162	REAR LAMPS BASE ASSY	48	32600494	REAR WHEEL ASSY, 16*4.00-10
17	44689070	OUTER TIRE. 14*3.00-8, BLACK, IA-1605	30	61200176	TURN LIGHT(L)	49	42500163	INNER REAR RIME ONLY, FOR 14*4.00-10
18	44600590	INNER TIRE. 3.00-8, BLACK	31	61200177	TURN LIGHT(R)	50	42500258	OUTTER RIM ONLY, FOR 14*4.00-10
19	53400023	BEARING ONLY, FORK 7 PCS/SET	32	61200178	REAR LIGHT	51	44600592	OUTER TIRE ONLY, 16*4.00-10, BLACK
20	63345605	WIRE, CONTROLLER	33	63400054	REAR LIGHT CABLE	52	44600593	INNER TIRE. 3.00-10, BLACK
21	64520006	CONTROLLER.	34	30519372	FRONT WINDSHIELD ONLY, SILVER WHITE	53	43601375	SHROUD, FENDER, LH
22	61230107	PG SD CONTROLLER PROTECTION KIT PG D51058	34	30561037	FRONT WINDSHIELD ONLY, ORANGE	54	43601377	SHROUD, FENDER, RH
23	30519553	FRONT SHROUD ASSY, SILVER WHITE	34	30550050	FRONT WINDSHIELD ONLY, SILVER DEEP GREEN	55	30500165	REAR SHROUD ASSY
23	30561981	FRONT SHROUD ASSY, ORANGE	34	30541079	FRONT WINDSHIELD ONLY, BRIGHT RED	56	30500346	SHROUD, MOTOR, BLACK
23	30550043	FRONT SHROUD ASSY, SILVER DEEP GREEN	35	30519400	FRONT WINDSHIELD FOLLOWING COVER ONLY, SILVER WHITE			



REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION	REF#	PART NO.	DESCRIPTION
57	32200388	CAPTAIN SEAT WITH HEADREST ONLY	69	43690171	THROTTLE LEVER ONLY	81	30541066	FRONT TILLER SHROUD ASSY, BRIGHT RED,S940
58	41700077	SLIDE MECH. ASSY., WITH RELEASE LEVER	70	61100069	POWER CONVERT BOARD	82	61200179	TURN LIGHT(L)
59	32100039	SEAT BASE WITH RELEASE LEVER ASSY	71	63680016	CABLE, POWER CONVERT	83	61200180	TURN LIGHT(R)
60	31000350	ARMREST ASSY., LH	72	63500036	HORN WIRE	84	30500365	REAR TILLER SHROUD ASSY
61	31000351	ARMREST ASSY., RH	73	61170007	HORN ONLY	85	33700256	CHARGER SOCKET SET
62	43700401	ARM PAD,BLACK	74	63180046	MAIN CABLE FOR TILLER	86	63700036	KEY SWITCH W/CONNECTING WIRE, S830
63	43700157	SAFETY BELT, 40" STD.	75	63400055	FRONT LIGHT CABLE	87	33700217	KEY SET
64	33500427	TILLER ASSY, SILVER	76	63742003	SENSOR, HALF SPEED	88	61138003	CABLE, PLUG.
64	33500428	TILLER ASSY,ORANGE	77	33500410	UPPER DASHBOARD ASSY	89	41700215	BRAKE CABLE ASSY
64	33500429	TILLER ASSY, SILVER DEEP GREEN	78	61020006	POTENTIOMETER., VR30K, COSMOS RV24	90	42700041	LEVER BRAKE, RH
64	33500430	TILLER ASSY, BRIGHT RED,S940	79	63500018	CONNECTING WIRE, SPEED POT	91	44689114	RUBBER BOOT
65	45700107	HANDLE BAR PAD	80	61108237	LCD PCB BOARD,S940	92	46781041	LEFT SIDE MIRROR,FO-297,2 TONE
66	41780053	ELEVATING CYLINDER	81	33500170	FRONT TILLER SHROUD ASSY., SILVER	93	46781042	RIGHT SIDE MIRROR,FO-297,2 TONE
67	63500021	CONNECTING WIRE (POT)	81	30561039	FRONT TILLER SHROUD ASSY,ORANGE	94	33700462	CHARGER,HP0240WB 24V 8A,WITH CORD
68	61010047	POTENTIOMETER	81	30550052	TILLER FRONT SHROUD ASSY, SILVER DEEP GREEN			