

Aluminium wheelchair L470-20/22/24

Features

Reinforced double cross-brace for extra strength with weight capacity of 182kg

Detachable full length armrests

Detachable swing away leg rests with oversize footplates, with optional elevating leg rests

24" x 1.3/8" puncture free foam filled tyres



Technical Specs

Overall Width 73.5cm(20") 78.5cm(22") 83.5cm(24")

Overall Length 123cm

Overall Height 91 cm

Seat Width 48cm(20") 56cm(22") 61cm(24")

Seat Depth 46cm

Seat height from ground 49cm

Seat Back Height 40cm

Footrest Length Adjustment 37cm – 50cm

Armrest Height 23cm

Weight Capacity 182kg

Total weight 15.5kg(20") 16.5kg(22") 17.5kg(24")