

# 4PT™

4 POLE TECHNOLOGY



## 745 PLUS

  
**AUSTRALIA**


[www.meritsaustralia.com](http://www.meritsaustralia.com)

# 745 PLUS

## Features

- Stylish eye catching design with attractive modern colour range
- This new design offers generous foot space & leg room, also with large battery capacity
- Loop handlebars offer comfortable finger or thumb throttle operation
- Adjustable high back contoured seat with simple operation of slide and swivel mechanisms to allow maximum user comfort
- Infinite adjustment of steering tiller allowing everyone to easily find a comfortable position
- High mount charger socket means easy access when charging
- 11" wheels which provide extra ground clearance
- Adjustable suspension to give maximum comfort
- 4 pole motor delivering more power
- Adjustable power saving LED headlight, indicators front and rear, rear tail lights and brake lights
- Built in speed reduction switch which is activated when turning offering exceptional safety.

## Dimensions cm/in



Overall Dimension A x B x C	132 x 61 x 102cm/52 x 24 x 40"
Seat Dimension : D x E x F	46 x 40 x 46cm/18 x 16 x 18"
Seat Height from Ground : G	61-66cm/24-26"
Seat Height from Deck : H	44-49cm/17-19"
Wheel Base : I	98cm/38"

## Weight kg/lb

base	battery/ea	seat	total with battery	capacity
51/112	16/35	15/33	98/217	150/330



## Technical Specs

Max Speed	10kph/6.25mph
Range	30-40km/20-26mi
Turn Radius	167cm/66"
Ground clearance	10cm/3.9"
Motor	DC 24V, 350W (Rating)/1200W (Max)
Controller	S-Drive 120A
Battery	12V/55 - 60A/H x 2pcs
Charger	4A off-board
Gradient	6°
Front Wheel	11" Pneumatic tire
Rear Wheel	11" Pneumatic tire
Brake	intelligent, regenerative, electromagnetic brakes

Australian Standards Compliant AS/NZS 3695.2:2013 Pending

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.