

Pull Apart Commode C324

Features

Pull Apart 3 in 1 commode can be used as a toilet safety frame, a raised toilet seat or outside the bathroom as a commode

Easy to clean aluminium finish

Durable plastic seat & lid, commode bucket, cover & splash shield

Armrest design provides added comfort



Technical Specs

Width between arms	46cm
Outside width	64cm
Seat width	36 cm
Seat Depth	38cm
Seat height	46-61cm
Overall height	66-81cm
Weight Capacity	136kg
Total weight	5kg