

Fende S946A

Features

- Stylish new modern design with a modesty panel gives this scooter a smart attractive look.
- Adjustable prestige car like suspension front & rear delivers a smooth ride over varied terrain.
- Comfortable leather feel contoured high back seat with slide, swivel, recline & adjustable armrests
- Unique battery storage design allows for greater floor space & a safe secure lockable storage area for personal items.
- New modern LCD dashboard display offers clearer more easy to read battery gauge, speed, time & temperature.
- LED lighting system delivers brighter illumination with less power consumption.
- Built in speed reduction switch which is activated when turning offering exceptional safety.
- With new easy to access 12 volt socket (for phone or GPS)



Lockable storage drawer underneath the seat.




Extra storage & Shopping bag hook.



Digital LCD Dashboard Display.

Dimensions cm/in



Overall Dimension : A x B x C 120 x 63 x 132cm/47.2 x 24.8 x 52"

Seat Dimension : D x E x F 46 x 41 x 46cm/18 x 16 x 18"

Seat Height from Ground : G 64-69cm/25-27"

Seat Height from Deck : H 41.5-46.5cm/16.3-18.3"

Wheel Base : I 88cm/34.6"

Weight kg/lb

Base	Battery/ea	Seat	Total with Battery	Capacity
70/154	6.5/14	18.5/40.5	114.5/252	136/300

Technical Specs

Max Speed	10kph/6.25mph
Range	30km/19mi
Turn Radius	143cm/56.3"
Ground clearance	10cm/3.9"
Motor	DC 24V, 295W (Rating)/1150W (Max)
Controller	Rhino 2 90A
Battery	12V/22A/H x 4pcs
Charger	5A off-board
Gradient	10°
Front Wheel	11"
Rear Wheel	11"
Brake	intelligent, regenerative, electromagnetic brakes

merits
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.