

# Rollator - W462

## Features

Lightweight aluminium frame

Durable 6" puncture free foam tyres

Easy to fold for storage & transportation

Padded seat with removable folding backrest



colour available - Green or Red

## Technical Specs

<b>Seat Width</b>	35cm
<b>Seat Depth</b>	33cm
<b>Seat Height</b>	53 cm
<b>Width at base</b>	64cm
<b>Depth at base</b>	61cm
<b>Height</b>	82-95cm
<b>Weight Capacity</b>	114kg
<b>Total weight</b>	6.4kg