



## Maverick 10




# Maverick 10

## ➤ Features

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors even over rough terrain
- Dual in-line motors enhanced efficiency, torque, range and performance
- Full front and rear suspension design offers a comfortable ride
- Low center of gravity design delivers a smooth and stable ride
- Swing back joystick mount is standard
- Height and width adjustable armrest, semi reclining backrest and height adjustable seat for various needs
- Available in a number of vibrant new colours – Red, Silver, Pearl Gold, Charcoal, NV Green, Kite Blue

## ➤ Dimensions cm / in



Overall Dimension : A x B x C	86 x 61 x 123cm / 34 x 24 x 48in
Seat Dimension : D x E x F	46 x 46 x 50cm / 18 x 18 x 20in
Seat Height from Ground : G	51–56cm / 20–22in
Seat Height from Deck : H	43–53cm / 17–21in

## ➤ Weight kg / lb

base	battery/ea	seat	total with battery	capacity
57 / 125	12 / 26	19 / 42	100kg / 220lb	136kg / 300lb



## ➤ Technical Specs.

Max Speed	8kph / 5mph
Range	28km / 18mi.
Turn Radius	50cm/ 20"
Ground clearance	5cm / 2"
Motor	DC24V, 160W(Rating) / 520W(Max)
Controller	PG VR2 40A
Battery	12V / 30-40AH x 2pcs
Charger	5A off-board
Gradient	8°
Caster	6" x 2" PU tire
Drive Wheel	10" (260x85) foam filled tire
Brake	intelligent, regenerative, electromagnetic brakes